Úvod

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Introduction

ENVY Therapy® HAIR Ultrasonic Therapy is a device for home use emitting ultrasonic waves and LED light effect on hair and scalp. You can choose from four mode settings (LED light colors) and three degrees of intensity.

Ultrasonic comb:

- Improves, stimulates and nourishes hair follicles.
- Reduces hair loss.
- Supports the growth and density of new hair.
- Reduces the symptoms of skin problems such as dandruff, excess grease and more.
- Improves the absorption of nutrients from products.
- Protects the scalp and maintains its proper pH.

- Read the instructions carefully before using the product.
- The comb is an electronic device for home use and is not suitable for diagnostic, therapeutic or other purposes.
- Thoroughly clean the device after use to prevent the accumulation of grease and other impurities.
- Fully charge before first use.

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INTENSIVE TREATMENT to reduce hair loss and support hair follicles

The red LED light (LED button - mode 2) regenerates hair follicles, increases their vitality and reduces hair loss.

The pink-red LED light (LED button - mode 3) supports hair growth and density, ensures freshness and protection of the scalp, fights local inflammations and strengthens the skin's immunity in the treated area.

Pulsing pink-red LED light (LED button - mode 4) improves nutrient absorption and nourishes hair follicles.

First 2 weeks - 1 ampoule twice a week with red LED light

3. - 6. week - 1 ampoule once a week with pink-red LED light

Use 1 ampoule every other week with a pulsing pink-red LED light as a maintenance treatment throughout the year.

SUPPORT TREATMENT to support hair growth and maintain the firmness and quality of the scalp

1 ampoule once a week for 6 weeks with a pink-red LED light.

Use 1 ampoule every other week with a pulsing pink-red LED light as a maintenance treatment throughout the year.

2.4. Charging

- Fully charge the comb before first use.
- Connect one end of the charging cable to the charging port on the comb and the other end to the adapter, plug it into an electrical socket. The product must not be used during charging.
- During charging, the charge level indicator will flash red.
- When charging is complete, the charge level indicator stays green.
- Charging takes about 2 and a half hours.
- The device lasts approximately 70 minutes when fully charged.
- Each therapy lasts 10 minutes, then the comb automatically shuts down.
- The number of uses may vary depending on the usage.
- Use the provided charging cable, or a charging cable that meets product specifications, and plug it into the adapter and electrical socket correctly to prevent malfunctions or fire caused by heat.

3. Warnings

- Read these warnings before using and storing the comb.
- Ultrasound treatment is not recommended for people with pacemakers or people suffering from epilepsy or severe heart disease.
- Do not store the device at too low or high temperatures, or in a humid environment. Store in the original box in a dry and cool place. Do not expose it to direct sunlight.
- Do not use the device while bathing or showering.

- For hygienic reasons, do not share this product with other people.
- Use the device only with ENVY Therapy® Hair Ampoules, as other products may not have the appropriate viscosity and may clog the internal teeth discharging solution, which is not covered by any product complaints.
- If you are using the device for the first time or if you have sensitive scalp, try it on a small skin area first and make sure your skin is in good condition before using the device.
- Check the device before use. If it is deformed, damaged, scratched, dirty or if it smells, do not continue to use it.
- Do not leave the product with the teeth of the comb facing downwards to prevent the solution from leaking out.
- If you have skin problems, consult your doctor about using the device.
- Immediately after using the device, a slight reddening of the skin may occur due to the intensity of the ultrasound (RP) and moisture of the skin.
- You may feel heat or tingling during use.
- If any of the following situations occur, consult your doctor before using the device:
- skin symptoms or possible skin side effects,
- symptoms that may affect fine motor skills,
- a recent surgery or cosmetic intervention in the head area,
- if you have taken medicine that may affect your skin, heat perception or fine motor skills.
- In case of difficulties, stop using the product immediately.
- Do not use on open wounds.
- The device should not be used by persons who are physically weak, have delayed reactions or mental disorders, unless they use it safely under supervision.
- The device must not be used by unattended children or people with disabilities. The device is not intended for playing.
- If the device becomes excessively hot during use, stop using it immediately to prevent burns.
- To reduce the risk of electric shock, short circuit, or fire, do not use a damaged charger or electrical socket.
- Do not use the device near flammable or explosive materials (e.g. alcohol, sprays, gasoline, etc.) as it may cause an explosion or fire.
- Unplug the device before cleaning and maintenance.
- Do not use the device for any purpose other than that specified in the instructions.
- The product is not waterproof and should not be immersed in water or rinsed with water.
- Use the device at an ambient temperature between 15 °C and 30 °C. Charge and store the device at an ambient temperature between 5 °C and 35 °C.
- Do not store the device in a humid environment.
- Arbitrary assembly, disassembly, repair and modification of the device can lead to failure, which will not be possible to repair or complain about. Do not attempt to repair the device and always contact the seller.
- Carry the comb in the original box to avoid damage.
- Do not recommend using the comb during pregnancy. Due to the complexity of the clinical trial and related serious ethical issues, the device has not been tested in pregnant women and therefore it is not recommended to be used during this period.

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